

RPA 及其应用简介

---新型核酸扩增技术



默瑞生物 市场部



- RPA背景



- RPA产品



- RPA应用

Recombinase Polymerase Amplification

重组酶 聚合酶 扩增技术

RPA是由多种酶和蛋白参与下，在恒温条件下实现核酸指数扩增的技术


- ◆ 高灵敏度，高特异性，反应快速
- ◆ 无需温度循环仪（themocycler）
- ◆ 非常适用于POCT, 疾病诊断, 流行病学监测, 食品安全检测, 动物病原物检测等

RPA---英国剑桥TwistDx公司专利技术

2006年

2010年

2014年



OPEN ACCESS Freely available online PLOS BIOLOGY

DNA Detection Using Recombination Proteins

Olaf Piepenburg¹, Colin H. Williams¹, Derek L. Stemple², Niall A. Armes^{1*}




1 KEM Scientific Ltd, Cambridge, United Kingdom, 2 Wellcome Trust Sanger Institute, Hinxton, Cambridge, United Kingdom




DNA amplification is essential to most nucleic acid testing strategies, but established techniques require sophisticated equipment or complex experimental procedures, and their uptake outside specialised laboratories has been limited. Our novel approach, recombinase polymerase amplification (RPA), couples isothermal recombinase-driven primer targeting of template material with strand-displacement DNA synthesis. It achieves exponential amplification with no need for pretreatment of sample DNA. Reactions are sensitive, specific, and rapid and operate at constant low temperature. We have also developed a probe-based detection system. Key aspects of the combined RPA amplification/detection process are illustrated by a test for the pathogen methicillin-resistant *Staphylococcus aureus*. The technology proves to be sensitive to fewer than ten copies of genomic DNA. Furthermore, products can be detected in a simple sandwich assay, thereby establishing an instrument-free DNA testing system. This unique combination of properties is a significant advance in the development of portable and widely accessible nucleic acid-based tests.

Citation: Piepenburg O, Williams CH, Stemple DL, Armes NA (2006) DNA detection using recombination proteins. *PLoS Biol* 4(7): e204. DOI: 10.1371/journal.pbio.094024



RPA—酶反应驱动的核酸扩增技术

 **Strand displacing polymerase**
 **Recombinase**
 **SSB**



核酸蛋白复合物合成




D环形成；单链固定




链置换；延伸

RPA---特异性高的探针检测系统



TwistAmp[®] exo system---Real time RPA



检测方法



TwistAmp nfo system---LFD RPA



RPA---低能耗的核酸扩增技术




Table 1 Characteristics of RPA and other isothermal amplification techniques

| Techniques | Templates | No. of Primers | No. of Enzymes | Pre-heating | Temperature (°C) | Time (min) | Products type |
|------------|--------------|------------------------|----------------|-------------|------------------|------------|---------------|
| RPA | DNA, RNA | 2 | 2 | No | 25 ~ 42 | 20 | DNA |
| LAMP | DNA | 4 ~ 6 | 1 | No | 60 ~ 65 | 60 ~ 90 | DNA |
| SDA | DNA | 4 | 2 | Yes | 37 | 120 | DNA |
| MDA | DNA | Random hexamer primers | 1 | No | 30 | 480 ~ 960 | DNA |
| RCA | Circular DNA | 1 | 2 | No | 37 | 60 | Circular DNA |
| HDA | DNA, RNA | 2 | 2 | No | 37, 60 ~ 65 | 120 | DNA |
| NASBA | DNA, RNA | 2 | 2 ~ 3 | Yes | 37 ~ 42 | 90 ~ 120 | RNA |
| SMART | DNA, RNA | 2 | 2 ~ 3 | Yes | 41 | 120 | DNA |

RPA---产品性能卓越，应用广泛

| | LAMP(环介导等温扩增)试剂盒 | RPA试剂盒 |
|------------------|------------------|----------------------|
| 主要组分 | 链置换DNA聚合酶 | 重组酶；单链结合蛋白；链置换DNA聚合酶 |
| 最适反应温度 | 60-65℃ | 37-40℃ |
| 推荐反应时间 | 40min左右 | 20min左右 |
| 引物数 | 2对 | 1对 |
| 保存形式 | 液体 | 冻干粉 |
| 检测方式 | 荧光目视；比浊法实时浊度仪 | 电泳；试纸条（肉眼观察）；荧光定量 |
| 扩增产物可否用于下游应用 | 否 | 是 |
| 可否Multiplex/多重反应 | 否 | 是 |

RPA---产品线齐全，最大限度满足客户需求

TwistAmp系列—供客户自行研发

Basic kit :
凝胶电泳检测

A



无需任何仪器耗材
凝胶电泳检测
适用样本：DNA,RNA

Exo kit :
荧光定量检测

B



需要荧光扩增仪
荧光定量检测
适用样本：DNA,RNA

Nfo kit :
侧流试纸条检测

C



需要胶体金试纸条
终点检测
适用样本：DNA

RPA---产品线齐全，最大限度满足客户需求



Twirla™ *mixing incubator*



Ball Dispenser MICRO



Bearing Balls

New!



T-8 Isothermal Device



T-16 Isothermal Device




通量：8 samples 净重：1.9kg
双通道：FAM; HEX
触屏操作，含混匀功能
仪器或者PC读取

通量：16 samples 净重：1.9kg
三通道：FAM; HEX ROX
触屏操作，含混匀功能
仪器或者PC读取

RPA---产品线齐全，最大限度满足客户需求



Milenia Hybridtech 1



侧流检测试纸条



Tough carry case



External battery

Real-time RPA检测结核分支杆菌 (MTB)

OPEN ACCESS Freely available online

PLOS ONE

Rapid Detection of *Mycobacterium tuberculosis* by Recombinase Polymerase Amplification

David S. Boyle^{1,3}, Ruth McNerney^{2,3}, Hwee Teng Low², Brandon Troy Leader¹, Ailyn C. Pérez-Osorio³, Jessica C. Meyer⁴, Denise M. O'Sullivan², David G. Brooks⁴, Olaf Piepenburg⁴, Matthew S. Forrest^{4*}

¹ Program for Appropriate Technology in Health, Seattle, WA, United States of America, ² Faculty of Infectious and Tropical Diseases, London School of Hygiene & Tropical Medicine, London, United Kingdom, ³ Washington State Department of Health, Public Health Laboratories, Shoreline, WA, United States of America, ⁴ TwistDx Limited, Cambridge, United Kingdom




Table 4. Tuberculosis detection by indirect smear microscopy vs RPA IS6110.

| | Culture positive | Culture negative |
|------------------------------------|------------------|------------------|
| Indirect smear microscopy positive | 34 | 5 |
| Indirect smear microscopy negative | 14 | 37 |
| RPA IS6110 positive | 42 | 2 |
| RPA IS6110 negative | 6 | 40 |

Testing pulmonary specimens (n = 90) by indirect smear microscopy and RPA IS6110 to detect tuberculosis, with comparison to liquid culture based test data. RPA IS6110 was more sensitive than indirect smear microscopy (87.5% (95% CI: 81.7, 93.2) vs 70.8% (95% CI: 62.91, 78.75)) and also more specific (95.4 (95% CI: 92.3, 98.1) vs 88% (95% CI: 83.6, 92.4)).

doi:10.1371/journal.pone.0103091.t004

Table 5. Tuberculosis detection by indirect smear microscopy vs RPA IS1081.

| | Culture positive | Culture negative |
|------------------------------------|------------------|------------------|
| Indirect smear microscopy positive | 31 | 4 |
| Indirect smear microscopy negative | 5 | 31 |
| RPA IS1081 positive | 32 | 0 |
| RPA IS1081 negative | 4 | 35 |

Testing pulmonary specimens (n = 71) by indirect smear microscopy and RPA IS1081 to detect tuberculosis, with comparison to liquid culture based test data. RPA IS1081 was more sensitive than indirect smear microscopy (91.4% (95% CI: 85.98, 9) vs 86.1% (95% CI: 78.1, 94.1)) and also more specific (100% vs 88.6% (95% CI: 80.8, 96.1)).

doi:10.1371/journal.pone.0103091.t005

RPA---病原体检测神器

RT-RPA检测禽流感病毒H7N9

Contents lists available at ScienceDirect


Journal of Clinical Virology



journal homepage: www.elsevier.com/locate/jcv

Diagnostics-in-a-Suitcase: Development of a portable and rapid assay for the detection of the emerging avian influenza A (H7N9) virus

Ahmed Abd El Wahed^{a,b,*}, Manfred Weidmann^c, Frank T. Hufert^d

^a Unit of Infection Models, German Primate Center, Kellnerweg 4, 37077 Goettingen, Germany
^b Department of Virology, Mansoura University, Algomhoria Street, 35516 Dakahlia, Egypt
^c Institute of Aquaculture, University of Stirling, FK9 4LA Stirling, Scotland, UK
^d Institute of Microbiology and Virology, Brandenburg Medical School Theodor-Fontane, Fehrbelliner Straße 38, 16816 Neuruppin, Brandenburg, Germany

 CrossMark




RPA---病原体检测神器

LFD-RPA检测沙眼衣原体

Sensitive and Rapid Detection of *Chlamydia trachomatis* by Recombinase Polymerase Amplification Directly from Urine Samples

Katrin Krölov,^{*} Jekaterina Frolova,^{*} Oana Tudoran,^{*†} Julia Suhorutsenko,^{*} Taavi Lehto,^{*} Hiljar Sibul,^{*} Imre Mäger,^{*} Made Laanpere,[‡] Indrek Tulp,^{§¶} and Ülo Langel^{*||}



RT-RPA检测埃博拉病毒

SCIENTIFIC REPORTS


OPEN Development and Evaluation of a Rapid and Sensitive EBOV-RPA Test for Rapid Diagnosis of Ebola Virus Disease

Received: 10 March 2016
Accepted: 11 May 2016
Published: 01 June 2016

Mingjuan Yang^{1,*}, Yuehua Ke^{1,2,*}, Xuesong Wang^{1,2,*}, Hang Ren^{1,*}, Wei Liu^{1,2,*}, Huijun Lu^{2,3,*}, Wenyi Zhang^{1,2}, Shiwei Liu⁴, Guohui Chang^{1,2}, Shuguang Tian^{1,2}, Lihua Wang^{2,5}, Liuyu Huang¹, Chao Liu^{1,2}, Ruifu Yang^{2,6,7} & Zeliang Chen^{1,2,8}




RPA---与其他技术结合，应用潜力巨大




Lab on a Chip

PAPER
View Article Online
View Journal | View Issue


Cite this: *Lab Chip*, 2016, 16, 2309

A centrifugal direct recombinase polymerase amplification (direct-RPA) microdevice for multiplex and real-time identification of food poisoning bacteria†

Multiple RPA+ 微流控




METHODOLOGY ARTICLE
Open Access

Optimizing illumina next-generation sequencing library preparation for extremely at-biased genomes

Samuel O Oyola^{1*}, Thomas D Otto¹, Yong Gu¹, Gareth Maslen¹, Magnus Manske¹, Susana Campino¹, Daniel J Turner², Bronwyn Machnis¹, Dominic P Kwiatkowski¹, Harold P Swerdlow¹ and Michael A Quall¹

RPA+NGS



RPA---核酸检测技术革命



总结：RPA是一种快速，灵敏，高效的核酸扩增系统，适用于多种应用场景，搭配其他分子检测技术，有望替代PCR的分子生物学垄断地位，掀起一场核酸检测技术革命。

RPA---It really works !

Thank you for your attention

更多
问题

RPA官方网站: <http://www.twistdx.co.uk/>

RPA中文网: <http://twistdx.com.cn/>

huyadong@moreybio.com

默瑞生物市场部 胡亚东 13382194023
